

PROFICIENCY SKILLS CHECKLIST 3: LARGE MEDICAL RESTRAINTS

KEY

POSITIONING	Roosevelt	IPS Car Seat	Wallaroo	Other (write in)
Weight/Height Guidelines	35-115 pounds 33.5-62 inches	20-102 pounds 36-60 inches	22-106 pounds Up to 56 inches	
Medical Conditions	Behavioral/developmental; obesity; neuromuscular; spina bifida; scoliosis	Behavioral/developmental; scoliosis; neuromuscular; spina bifida	Neuromuscular; behavior; obesity; hip casts	
Position Dolls in Restraint	Head support; recline bar required under 65 pounds; at or above shoulders	Set up seat using plastic adjusters; do not unthread, just pass through slots; head support; at or above shoulders	Harness pads at all times; at or above shoulders	
Adjust Harness to Fit Properly	European harness routing; 2 nd crotch strap can be used after 65 pounds	Fine adjusters on front of seat	A-lock adjuster	
Adjust Harness Clip Properly	Mid-chest; can get anti-escape chest clip	Mid-chest	Mid-chest	
INSTALLATION	Roosevelt	IPS Car Seat	Wallaroo	Other (write in)
Tether Recommendations/Limits	EZ-Tether preferred; tethers not attached; 2 tethers over 80 pounds; can order HD anchors; call about tethering to back seat belt	Must tether; 2 tethers over 80 pounds; can order HD anchors; call about tethering to back seat belt	Must tether; 2 tethers over 40 pounds; call about tethering to back seat belt	
Lower Anchor Recommendations/Limits	N/A	N/A	Can be used 22-40 pounds	
Installed This Restraint	Pull harness tight; long belt path with EZ-Tether or top tether strap; short belt path with top tether strap(s)	Seat belt over metal frame; seat cannot extend more than 1 inch over vehicle seat	LATCH: 22-40 pounds Lap belt: 22-40 pounds (1 tether); 40- 106 pounds (2 tethers) Lap/shoulder belt: 22-50 pounds (1 tether); 50-106 pounds (2 tethers)	
INSTRUCTOR SIGN-OFF				