

PROFICIENCY SKILLS CHECKLIST 2: LARGE MEDICAL RESTRAINTS

NAME:					DATE:
POSITIONING	Roosevelt	Spirit/Spirit Plus	Wallaroo (now discontinued)	Other (write in)	Other (write in)
Weight/Height guidelines	35-115 pounds 33.5-62 inches	25-130 pounds Up to 66 Inches	22-106 pounds Up to 56 inches		
Medical Conditions	Behavioral/developmental; obesity; neuromuscular; spina bifida; scoliosis	Behavioral/developmental; scoliosis; neuromuscular; spina bifida	Neuromuscular; behavior; obesity; hip casts		
Position Dolls in restraint	Head support; recline bar required under 65 pounds; at or above shoulders	Set up seat using adjusters; do not unthread, just pass through slots; at or above shoulders	Harness pads at all times; at or above shoulders		
Adjust harness to fit properly	European harness routing; 2 nd crotch strap can be used after 65 pounds	Fine adjusters on front of seat, Spirit Plus swing away laterals and hip guides need adjusting	A-lock adjuster		
Adjust harness clip properly	Mid-chest; can get anti- escape chest clip	Mid-chest	Mid-chest		
INSTALLATION	Roosevelt	Spirit/Spirit Plus	Wallaroo (now discontinued)	Other (write in)	Other (write in)
Tether recommendations/limits Type/Preinstalled/HD	EZ-Tether preferred; tethers not attached; 2 tethers over 80 pounds; can order HD anchors; call about tethering to back seat belt	Must tether; 2 tethers over 80 pounds; can order HD anchors; call about tethering to back seat belt	Must tether; 2 tethers over 40 pounds; call about tethering to back seat belt		
LATCH recommendations/limits	N/A	N/A	Can be used 22-40 pounds		
Installed this restraint	Pull harness tight; long belt path with EZ-Tether or top tether strap; short belt path with top tether strap(s)	Seat belt with long belt path; Must tether; 2 tethers over 80 pounds	LATCH: 22-40 pounds Lap belt: 22-40 pounds (1 tether); 40-106 pounds (2 tethers) Lap/shoulder belt: 22-50 pounds (1 tether); 50-106 pounds (2 tethers)		